



Free Standing

Hydrotherapy Pools

www.HYDROPHYSIO.com



Gently Maintaining fitness and muscle tone

Hydrotherapy assists in post-operative rehabilitation, pain management, weight loss and psychological wellbeing. Low impact exercise is gentle on ligaments and joints, helping to keep patients as comfortable and pain free as possible.

This free standing pool system is recommended for canine swimming. The liner of the pool is manufactured from high tenacity polyester fabric, coated with PVC on both sides. The material resistance is greater than 220 or 280kg/5cm (depending on the model).

The anti-UVA and fungicide properties of the liner material ensures durability and creates a clean, safe environment. The pool frames are anodised aluminium, which is robust, rigid and lightweight. When assembled with the cast aluminium joints, the pools will stand upright without water. The pools are supplied with comprehensive assembly instructions and we can recommend installation teams if you require assistance during setup.



Pool Details and options




Pool Details:

- Kit format pools for easy installation and relocation
- No groundworks required
- Reinforced pool liner developed for canine use
- A low cost hydrotherapy option
- Water treatment systems, recirculation including a skimmer

Options:

- Five pool sizes (see rear page)
- Non-slip external and internal patient access ramps
- Non-slip internal step unit
- Therapist 'up and over' ladder
- Resistance jets
- Conventional electric or air source heat pump for heating the water
- Pool cover

Models and sizes



MAX - 5

56m³

1165 x 490cm

MAX - 4

50m³

1040 x 490cm

MAX - 3

32m³

915 x 365cm

MAX - 2

23m³

665 x 365cm

MAX - 1

18m³

540 x 365cm

All models are 120cm in depth and measurements given are external

See Our HYDRO PHYSIO Treadmills

The ideal combination for hydrotherapy suites. Systems can be installed poolside, utilising the pool water as a source for the treadmill unit.



FM 25162

MD 516639

OHS 672496